

Summary

Production Name	CYP2A7 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CYP2A7
Alternative Names	CYP2A7; Cytochrome P450 2A7; CYP11A7; Cytochrome P450 11A4
Gene ID	1549.0
SwissProt ID	P20853.Synthesized peptide derived from the C-terminal region of human CYP2A7.

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:40000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	56kDa

Background

Product Name: CYP2A7 Rabbit Polyclonal Antibody
Catalog #: APRab09649




This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum; its substrate has not yet been determined. This gene, which produces two transcript variants, is part of a large cluster of cytochrome P450 genes from the CYP2A, CYP2B and CYP2F subfamilies on chromosome 19q. [provided by RefSeq, Jul 2008], catalytic activity: $\text{RH} + \text{reduced flavoprotein} + \text{O}_2 = \text{ROH} + \text{oxidized flavoprotein} + \text{H}_2\text{O}$, caution: The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., cofactor: Heme group., function: Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acids, and xenobiotics., induction: P450 can be induced to high levels in liver and other tissues by various foreign compounds, including drugs, pesticides, and carcinogens., similarity: Belongs to the cytochrome P450 family.,


Research Area

Caffeine metabolism; Retinol metabolism; Drug metabolism; Drug metabolism;

Image Data



Western Blot analysis of various cells using CYP2A7 Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of 22RV1 cells using CYP2A7 Polyclonal Antibody diluted at 1: 2000

Product Name: CYP2A7 Rabbit Polyclonal Antibody
Catalog #: APRab09649



Note

For research use only.