
Product Name: DLEU7 Rabbit Polyclonal Antibody**Catalog #: APRab10014**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:50-1:200,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000

Molecular Weight

Antigen Information

Gene Name	DLEU7 LEU7
Alternative Names	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
Gene ID	220107.0
SwissProt ID	Q6UYE1
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human DLEU7. AA range:51-100


Background

alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue

specificity:Ubiquitous.,alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue specificity:Ubiquitous.,

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200