

Product Name: DQX1 Rabbit Polyclonal Antibody**Catalog #: APRab10148**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Monkey
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:10000-1:20000
Molecular Weight	79kDa

Antigen Information

Gene Name	DQX1
Alternative Names	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1
Gene ID	165545.0
SwissProt ID	Q8TE96
Immunogen	The antiserum was produced against synthesized peptide derived from human DQX1. AA range:551-600


Background

similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,similarity:Contains 1


helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,

Research Area

Image Data



Western blot analysis of lysates from COS7 cells, using DQX1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using DQX1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .