
Product Name: ECP Rabbit Polyclonal Antibody**Catalog #: APRab10285**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	18kDa

Antigen Information

Gene Name	RNASE3
Alternative Names	RNASE3; ECP; RNS3; Eosinophil cationic protein; ECP; Ribonuclease 3; RNase 3
Gene ID	6037.0
SwissProt ID	P12724
Immunogen	The antiserum was produced against synthesized peptide derived from human ECP. AA range:102-151


Background

Human,Mouse,Rat,Rabbit,Fish,Chicken,Guinea pig,Sheep,Bovine

Research Area

Asthma;

Image Data



Western blot analysis of lysate from 293 cells, using ECP antibody.

Western Blot analysis of various cells using ECP Polyclonal Antibody