

---

**Product Name: ECRG4 Rabbit Polyclonal Antibody****Catalog #: APRab10286**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	IHC,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

**Dilution Ratio** IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000

**Molecular Weight**

**Antigen Information**

<b>Gene Name</b>	C2orf40
<b>Alternative Names</b>	C2orf40; ECRG4; Augurin; Esophageal cancer-related gene 4 protein
<b>Gene ID</b>	84417.0
<b>SwissProt ID</b>	Q9H1Z8
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human ECRG4. AA range:41-90


**Background**

function:Probable hormone.,similarity:Belongs to the augurin family.,function:Probable hormone.,similarity:Belongs to the


augurin family.,

## Research Area

### Image Data



Immunofluorescence analysis of A549 cells, using ECRG4 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using ECRG4 Antibody. The picture on the right is blocked with the synthesized peptide.