
Product Name: FoxB2 Rabbit Polyclonal Antibody**Catalog #: APRab11075**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:20000-1:40000
Molecular Weight	46kDa

Antigen Information

Gene Name	FOXB2
Alternative Names	FOXB2; Forkhead box protein B2
Gene ID	442425.0
SwissProt ID	Q5VYV0
Immunogen	The antiserum was produced against synthesized peptide derived from human FOXB2. AA range:251-300


Background

function:Transcription factor, similarity:Contains 1 fork-head DNA-binding domain, function:Transcription


factor.,similarity:Contains 1 fork-head DNA-binding domain.,

Research Area

Image Data



Western blot analysis of lysates from HeLa, LOVO, and A549 cells, using FOXB2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using FoxB2 Polyclonal Antibody.