

Product Name: GADL1 Rabbit Polyclonal Antibody
Catalog #: AP Rab11252



Summary

Production Name	GADL1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC-P, WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GADL1
Alternative Names	Glutamate decarboxylase-like protein 1 (EC 4.1.1.-)
Gene ID	339896.0
SwissProt ID	Q6ZQY3.

Application

Dilution Ratio	WB 1:500-2000, IHC-P 1:50-200
Molecular Weight	57kDa

Background


May catalyze the decarboxylation of L-aspartate, 3-sulfinyl-L-alanine (cysteine sulfinic acid), and L-cysteate to beta-alanine, hypotaurine and taurine, respectively. Does not exhibit any decarboxylation activity toward glutamate.

Product Name: GADL1 Rabbit Polyclonal Antibody
Catalog #: APRab11252



Research Area

Image Data



Western Blot analysis of Hela cell lysis, using primary antibody at 1:1000 dilution.

Note

For research use only.