

Summary

Production Name	GPR110 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GPR110
Alternative Names	GPR110; PGR19; Probable G-protein coupled receptor 110; G-protein coupled receptor KPG_012; G-protein coupled receptor PGR19
Gene ID	266977.0
SwissProt ID	Q5T601.The antiserum was produced against synthesized peptide derived from human GPR110. AA range:831-880

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:40000.Not yet tested in other applications.
-----------------------	---

Product Name: GPR110 Rabbit Polyclonal Antibody
Catalog #: AP Rab11625



Molecular Weight 110kDa


Background

function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7


subfamily., similarity: Contains 1 GPS domain., similarity: Contains 1 SEA domain., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily., similarity: Contains 1 GPS domain., similarity: Contains 1 SEA domain.,

Research Area

Image Data




Immunofluorescence analysis of A549 cells, using GPR110 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from Jurkat, HUVEC, and COLO cells, using GPR110 Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: GPR110 Rabbit Polyclonal Antibody
Catalog #: APRab11625



Western Blot analysis of various cells using GPR110 Polyclonal Antibody



Western Blot analysis of 3T3 cells using GPR110 Polyclonal Antibody

Note

For research use only.