
Product Name: GRIN3 Rabbit Polyclonal Antibody**Catalog #: APRab11768**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000

Molecular Weight

Antigen Information


Gene Name	GPRIN3
Alternative Names	GPRIN3; KIAA2027; G protein-regulated inducer of neurite outgrowth 3; GRIN3
Gene ID	285513.0
SwissProt ID	Q6ZVF9
Immunogen	The antiserum was produced against synthesized peptide derived from human GPRIN3. AA range:131-180

Background

function:May be involved in neurite outgrowth.,function:May be involved in neurite outgrowth.,

Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using GPRIN3 Antibody. The picture on the right is blocked with the synthesized peptide.