

Product Name: GUSBP1 Rabbit Polyclonal Antibody**Catalog #: APRab11861**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:20000-1:40000
Molecular Weight	18kDa

Antigen Information


Gene Name	GUSBP1
Alternative Names	GUSBP1; SMA3; Putative beta-glucuronidase-like protein SMA3
Gene ID	153561.0
SwissProt ID	Q15486
Immunogen	The antiserum was produced against synthesized peptide derived from human GUSBL1. AA range:51-100

Background


Binds GTP, GDP and GMP.skeletal system development,

Research Area

Image Data



Western blot analysis of the lysates from HT-29 cells using GUSBL1 antibody.



Western Blot analysis of K562 cells using GUSBP1 Polyclonal Antibody