

Product Name: IFI27 Rabbit Polyclonal Antibody**Catalog #: APRab12373**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	13kDa

Antigen Information


Gene Name	IFI27
Alternative Names	
Gene ID	3429.0
SwissProt ID	P40305
Immunogen	Synthesized peptide derived from human protein . at AA range: 1-80

Background

induction:By alpha interferons.,similarity:Belongs to the IFI6/IFI27 family.,induction:By alpha interferons.,similarity:Belongs to the IFI6/IFI27 family.,

Research Area

Image Data



Western blot analysis of lysates from HCT116 cells, IFI27 Rabbit Polyclonal Antibody was diluted at 1:1000, 4°over night