
Product Name: IgD Rabbit Polyclonal Antibody**Catalog #: APRab12419**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:50-1:200,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000

Molecular Weight

Antigen Information

Gene Name	IGHD
Alternative Names	Ig delta chain C region
Gene ID	
SwissProt ID	P01880
Immunogen	Synthetic peptide from human protein at AA range: 111-160


Background

similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,similarity:Contains 3 Ig-like (immunoglobulin-like) domains.,

Research Area

Immunology

Image Data



Immunohistochemical analysis of paraffin-embedded human-tonsil, antibody was diluted at 1:200