
Product Name: IGFL2 Rabbit Polyclonal Antibody**Catalog #: APRab12439**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	13kDa

Antigen Information


Gene Name	IGFL2 UNQ645/PRO1275
Alternative Names	
Gene ID	147920.0
SwissProt ID	Q6UWQ7
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50

Background

IGF like family member 2(IGFL2) Homo sapiens IGFL2 belongs to the insulin-like growth factor (IGF; see MIM 147440) family of signaling molecules that play critical roles in cellular energy metabolism and in growth and development, especially prenatal growth (Emtage et al., 2006 [PubMed 16890402]).[supplied by OMIM, Mar 2008],similarity:Belongs to the IGFL family.,

Research Area

Image Data



Western blot analysis of lysates from SH-SY5Y cells, IGFL2 Rabbit Polyclonal Antibody was diluted at 1:1000, 4°over night