
Product Name: IL-17C Rabbit Polyclonal Antibody**Catalog #: APRab12506**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:50-1:200,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000

Molecular Weight

Antigen Information

Gene Name	IL17C UNQ561/PRO1122
Alternative Names	Interleukin-17C (IL-17C;Cytokine CX2)
Gene ID	27189.0
SwissProt ID	Q9P0M4
Immunogen	Synthetic peptide from human protein at AA range: 81-130


Background

The protein encoded by this gene is a T cell-derived cytokine that shares the sequence similarity with IL17. This cytokine was reported to stimulate the release of tumor necrosis factor alpha and interleukin 1 beta from a monocytic cell line. The


expression of this cytokine was found to be restricted to activated T cells. [provided by RefSeq, Jul 2008],function:Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.,online information:Interleukin-17 entry,similarity:Belongs to the IL-17 family.,

Research Area


Image Data



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-colon, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-pancreas, antibody was diluted at 1:200