
Product Name: ISK7 Rabbit Polyclonal Antibody**Catalog #: APRab12771**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio WB 1:500-1:2000

Molecular Weight


Antigen Information

Gene Name	SPINK7
Alternative Names	Serine protease inhibitor Kazal-type 7 (Esophagus cancer-related gene 2 protein) (ECRG-2)
Gene ID	84651.0
SwissProt ID	P58062
Immunogen	Synthesized peptide derived from human ISK7 AA range: 1-50

Background

Research Area

Image Data



Western Blot analysis of hela 293T using primary antibody at 1:1000 dilution 4°C, overnight. Secondary antibody was diluted at 1:10000 25°C, 1.5hours