

Product Name: LDLRAD3 Rabbit Polyclonal Antibody
Catalog #: AP Rab13271



Summary

Production Name	LDLRAD3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	LDLRAD3
Alternative Names	LDLRAD3; Low-density lipoprotein receptor class A domain-containing protein 3
Gene ID	143458.0
SwissProt ID	Q86YD5. The antiserum was produced against synthesized peptide derived from human LDLRAD3. AA range:37-86

Application

Dilution Ratio	WB 1:500-2000, IHC-P 1:50-300
Molecular Weight	37kDa

Background


Product Name: LDLRAD3 Rabbit Polyclonal Antibody
Catalog #: APRab13271



similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,similarity:Belongs to the LDLR family.,similarity:Contains 3 LDL-receptor class A domains.,

Research Area

Image Data



Immunofluorescence analysis of A549 cells, using LDLRAD3 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.