

Product Name: mAChR M4 Rabbit Polyclonal Antibody
Catalog #: APRab13545



Summary

Production Name	mAChR M4 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CHRM4
Alternative Names	CHRM4; Muscarinic acetylcholine receptor M4
Gene ID	1132.0
SwissProt ID	P08173.The antiserum was produced against synthesized peptide derived from human CHRM4. AA range:236-285

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:10000.Not yet tested in other applications.
Molecular Weight	53kDa

Product Name: mAChR M4 Rabbit Polyclonal Antibody
Catalog #: APRab13545




Background

The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine and includes cellular responses such as adenylate cyclase inhibition, phosphoinositide degradation, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The clinical implications of this receptor are unknown; however, mouse studies link its function to adenylyl cyclase inhibition. [provided by RefSeq, Jul 2008],function:The muscarinic acetylcholine receptor mediates various cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels through the action of G proteins. Primary transducing effect is inhibition of adenylate cyclase.,similarity:Belongs to the G-protein coupled receptor 1 family.,


Research Area

Neuroactive ligand-receptor interaction;Regulates Actin and Cytoskeleton;

Image Data



Immunofluorescence analysis of LOVO cells, using CHRM4 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using mAChR M4 Polyclonal Antibody

Product Name: mAChR M4 Rabbit Polyclonal Antibody
Catalog #: APRab13545



Note

For research use only.