

Product Name: Mammaglobin A Rabbit Polyclonal Antibody
Catalog #: APRab13614

Summary

Production Name	Mammaglobin A Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	SCGB2A2
Alternative Names	SCGB2A2; MGB1; UGB2; Mammaglobin-A; Mammaglobin-1; Secretoglobin family 2A member 2
Gene ID	4250.0
SwissProt ID	Q13296.The antiserum was produced against synthesized peptide derived from human Mammaglobin. AA range:41-90

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:5000.Not
-----------------------	---

Product Name: Mammaglobin A Rabbit Polyclonal Antibody
Catalog #: APRab13614

yet tested in other applications.


Molecular Weight 20kDa

Background


similarity:Belongs to the secretoglobin family. Lipophilin subfamily.,tissue specificity:Mammary gland specific. Over-expressed in breast cancer.,similarity:Belongs to the secretoglobin family. Lipophilin subfamily.,tissue specificity:Mammary gland specific. Over-expressed in breast cancer.,

Research Area

Image Data




Immunofluorescence analysis of HepG2 cells, using Mammaglobin Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using Mammaglobin Antibody. The picture on the right is blocked with the synthesized peptide.

**Product Name: Mammaglobin A Rabbit Polyclonal
Antibody**
Catalog #: APRab13614



Western blot analysis of lysates from HepG2 cells, using Mammaglobin Antibody. The lane on the right is blocked with the synthesized peptide.

Note

For research use only.