

Product Name: MEGF6 Rabbit Polyclonal Antibody**Catalog #: APRab13788**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	169kDa

Antigen Information


Gene Name	MEGF6 EGFL3 KIAA0815
Alternative Names	
Gene ID	1953.0
SwissProt ID	O75095
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120

Background

similarity:Contains 25 EGF-like domains.,similarity:Contains 25 EGF-like domains.,

Research Area

Image Data



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution.
Secondary antibody was diluted at 1:10000