

Summary

Production Name	MRGF Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	MRGPRF MRGPRF; GPR140; GPR168; MRGF; PSEC0142; Mas-related G-protein coupled receptor member F; Mas-related gene F protein; G-protein coupled receptor 140; G-protein coupled receptor 168
Alternative Names	
Gene ID	116535.0
SwissProt ID	Q96AM1.The antiserum was produced against synthesized peptide derived from human MRGRF. AA range:271-320

Application

Dilution Ratio	WB 1:500-2000, IHC-P 1:50-300
Molecular Weight	40kDa

Product Name: MRGF Rabbit Polyclonal Antibody
Catalog #: AP Rab14087




Background

function: Orphan receptor. May bind to a neuropeptide and may regulate nociceptor function and/or development, including the sensation or modulation of pain., similarity: Belongs to the G-protein coupled receptor 1


family., similarity: Belongs to the G-protein coupled receptor 1 family. Mas subfamily., function: Orphan receptor. May bind to a neuropeptide and may regulate nociceptor function and/or development, including the sensation or modulation of pain., similarity: Belongs to the G-protein coupled receptor 1 family., similarity: Belongs to the G-protein coupled receptor 1 family. Mas subfamily.,

Research Area

Image Data



Western blot analysis of lysates from HUVEC cells, using MRGRF Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using MRGF Polyclonal Antibody

Product Name: MRGF Rabbit Polyclonal Antibody
Catalog #: APRab14087



Note

For research use only.