

**Product Name: Olfactory receptor 10D4 Rabbit Polyclonal Antibody**  
**Catalog #: APRab15133**

---

## Summary

|                        |                                                    |
|------------------------|----------------------------------------------------|
| <b>Production Name</b> | Olfactory receptor 10D4 Rabbit Polyclonal Antibody |
| <b>Description</b>     | Rabbit Polyclonal Antibody                         |
| <b>Host</b>            | Rabbit                                             |
| <b>Application</b>     | WB,IF-P,IF-F,ICC/IF,ELISA                          |
| <b>Reactivity</b>      | Human,Rat,Mouse                                    |

## Performance

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                                      |
| <b>Modification</b> | Unmodified                                                                                        |
| <b>Isotype</b>      | IgG                                                                                               |
| <b>Clonality</b>    | Polyclonal                                                                                        |
| <b>Form</b>         | Liquid                                                                                            |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| <b>Buffer</b>       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| <b>Purification</b> | Affinity purification                                                                             |

## Immunogen

|                          |                                                                                                           |
|--------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | OR10D4P                                                                                                   |
| <b>Alternative Names</b> | OR10D4P; OR10D4; Putative olfactory receptor 10D4                                                         |
| <b>Gene ID</b>           | 390266.0                                                                                                  |
| <b>SwissProt ID</b>      | Q8NGN7.The antiserum was produced against synthesized peptide derived from human OR10D4. AA range:161-210 |

## Application

|                       |                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------|
| <b>Dilution Ratio</b> | WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:10000.Not yet tested in other applications. |
|-----------------------|-----------------------------------------------------------------------------------------------------|

**Product Name: Olfactory receptor 10D4 Rabbit  
Polyclonal Antibody  
Catalog #: APRab15133**

---


**Molecular Weight** 38kDa

## Background


Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],caution:Could be the product of a pseudogene.,function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,

## Research Area

## Image Data




Immunofluorescence analysis of LOVO cells, using OR10D4 Antibody. The picture on the right is blocked with the synthesized peptide.




**Product Name: Olfactory receptor 10D4 Rabbit  
Polyclonal Antibody  
Catalog #: AP Rab15133**



Western blot analysis of lysates from COLO cells, using OR10D4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using OR10D4 antibody.



Western Blot analysis of various cells using Olfactory receptor 10D4 Polyclonal Antibody

## Note

For research use only.