

Product Name: Olfactory receptor 52E5 Rabbit Polyclonal Antibody
Catalog #: APRab15264

Summary

Production Name	Olfactory receptor 52E5 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR52E5
Alternative Names	OR52E5; Olfactory receptor 52E5
Gene ID	390082.0
SwissProt ID	Q8NH55.The antiserum was produced against synthesized peptide derived from human OR52E5. AA range:199-248

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:20000.Not yet tested in other applications.
-----------------------	---

Product Name: Olfactory receptor 52E5 Rabbit Polyclonal Antibody
Catalog #: APRab15264


Molecular Weight 36kDa

Background


Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor „similarity:Belongs to the G-protein coupled receptor 1 family,„

Research Area

Image Data



Western blot analysis of lysates from COLO, MCF-7, and HeLa cells, using OR52E5 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using OR52E5 antibody.

**Product Name: Olfactory receptor 52E5 Rabbit
Polyclonal Antibody
Catalog #: APRab15264**



Western Blot analysis of various cells using Olfactory receptor 52E5 Polyclonal Antibody

Note

For research use only.