

Product Name: Olfactory receptor 5H1 Rabbit Polyclonal Antibody**Catalog #: AP Rab15287**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB, ICC/IF, ELISA
Reactivity	Human, Rat, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
Molecular Weight	38kDa

Antigen Information

Gene Name	OR5H1
Alternative Names	OR5H1; Olfactory receptor 5H1; HTPCRX14
Gene ID	26341.0
SwissProt ID	A6NKK0
Immunogen	The antiserum was produced against synthesized peptide derived from human OR5H1. AA range:262-311


Background

olfactory receptor family 5 subfamily H member 1 (OR5H1) Homo sapiens Olfactory receptors interact with odorant molecules


in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor „similarity:Belongs to the G-protein coupled receptor 1 family.,

Research Area


Image Data



Immunofluorescence analysis of A549 cells, using OR5H1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using OR5H1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using OR5H1 antibody.