

Product Name: Olfactory receptor 5K3 Rabbit Polyclonal Antibody
Catalog #: APRab15292

Summary

Production Name	Olfactory receptor 5K3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR5K3
Alternative Names	OR5K3; Olfactory receptor 5K3
Gene ID	403277.0
SwissProt ID	A6NET4.The antiserum was produced against synthesized peptide derived from human OR5K3. AA range:232-281

Application

Dilution Ratio	WB 1:500-1:2000, IF-P/IF-F/ICC/IF 1:200-1:1000, ELISA 1:5000.Not yet tested in other applications.
-----------------------	--

Product Name: Olfactory receptor 5K3 Rabbit Polyclonal Antibody
Catalog #: AP Rab15292


Molecular Weight 37kDa

Background


olfactory receptor family 5 subfamily K member 3(OR5K3) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor „similarity:Belongs to the G-protein coupled receptor 1 family.,

Research Area

Image Data



Western blot analysis of lysates from K562 cells, using OR5K3 Antibody. The lane on the right is blocked with the synthesized peptide.





**Product Name: Olfactory receptor 5K3 Rabbit Polyclonal
Antibody**
Catalog #: APRab15292

Western blot analysis of the lysates from HepG2 cells using OR5K3 antibody.

Note

For research use only.