
Product Name: OSR2 Rabbit Polyclonal Antibody**Catalog #: APRab15522**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:20000-1:40000
Molecular Weight	35kDa

Antigen Information

Gene Name	OSR2
Alternative Names	OSR2; Protein odd-skipped-related 2
Gene ID	116039.0
SwissProt ID	Q8N2R0
Immunogen	Synthesized peptide derived from OSR2 . at AA range: 60-140

Background


odd-skipped related transcription factor 2(OSR2) Homo sapiens OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).[supplied by OMIM, Mar 2008],similarity:Belongs

to the Odd C2H2-type zinc-finger protein family.,similarity:Contains 5 C2H2-type zinc fingers.,

Research Area

Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger; Stem Cells; Mesenchymal Stem Cells; Osteogenesis

Image Data



Western Blot analysis of JK cells using OSR2 Polyclonal Antibody.