

Product Name: Pan myristoylation Rabbit Polyclonal Antibody**Catalog #: APRab15722**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Species independent
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:10000-1:20000
Molecular Weight	130 80 35kDa


Antigen Information

Gene Name	
Alternative Names	
Gene ID	
SwissProt ID	
Immunogen	Synthetic myristoylated peptide

Background

Research Area

Image Data



Western Blot analysis of KB cells using Antibody diluted at 500. Secondary antibody was diluted at 1:20000