

Product Name: Photomedin-2 Rabbit Polyclonal Antibody**Catalog #: APRab16084**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
Molecular Weight	84kDa

Antigen Information

Gene Name	OLFML2B
Alternative Names	OLFML2B; Olfactomedin-like protein 2B; Photomedin-2
Gene ID	25903.0
SwissProt ID	Q68BL8
Immunogen	The antiserum was produced against synthesized peptide derived from human OLFML2B. AA range:171-220


Background

PTM:O-glycosylated and N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin


and chondroitin sulfate E.,PTM:O-glycosylated and N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin and chondroitin sulfate E.,

Research Area

Image Data



Western blot analysis of lysates from HUVEC cells, using OLFML2B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Photomedin-2 Polyclonal Antibody