

Product Name: Pmp24 Rabbit Polyclonal Antibody
Catalog #: APRab16309



Summary

Production Name	Pmp24 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	PXMP4
Alternative Names	PXMP4; PMP24; Peroxisomal membrane protein 4; 24 kDa peroxisomal intrinsic membrane protein
Gene ID	11264.0
SwissProt ID	Q9Y6I8.The antiserum was produced against synthesized peptide derived from human PXMP4. AA range:15-64

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:20000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	24kDa

Product Name: Pmp24 Rabbit Polyclonal Antibody
Catalog #: APRab16309




Background


similarity:Belongs to the peroxisomal membrane protein PXMP2/4 family.,subunit:Interacts with PEX19.,similarity:Belongs to the peroxisomal membrane protein PXMP2/4 family.,subunit:Interacts with PEX19.,

Research Area

Image Data



Western blot analysis of lysates from COS7 and COLO cells, using PXMP4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of COS-7 cells using Pmp24 Polyclonal Antibody

Note

For research use only.