

**Product Name: PPP1R14D Rabbit Polyclonal Antibody****Catalog #: APRab16424**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC,ICC/IF,ELISA
<b>Reactivity</b>	Human,Rat,Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
<b>Molecular Weight</b>	24kDa

**Antigen Information**


<b>Gene Name</b>	PPP1R14D
<b>Alternative Names</b>	PPP1R14D; GBPI; Protein phosphatase 1 regulatory subunit 14D; Gastrointestinal and brain-specific PP1-inhibitory protein 1; GBPI-1
<b>Gene ID</b>	54866.0
<b>SwissProt ID</b>	Q9NXH3
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human PPP1R14D. AA range:66-115

**Background**


Protein phosphatase-1 (PP1; see MIM 176875) is a major cellular phosphatase that reverses serine/threonine protein phosphorylation. PPP1R14D is a PP1 inhibitor that itself is regulated by phosphorylation (Liu et al., 2004 [PubMed 12974676]). [supplied by OMIM, Feb 2010],function:Inhibitor of PPP1CA. Has inhibitory activity only when phosphorylated, creating a molecular switch for regulating the phosphorylation status of PPP1CA substrates and smooth muscle contraction.,PTM:Phosphorylated on several residues.,similarity:Belongs to the PP1 inhibitor family.,tissue specificity:Detected in colon, intestine, kidney and brain cortex.,

## Research Area

## Image Data



Western blot analysis of lysates from HeLa, HepG2, and COLO cells, using PPP1R14D Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using PPP1R14D Polyclonal Antibody