
Product Name: Rab 7L1 Rabbit Polyclonal Antibody**Catalog #: APRab16763**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
Molecular Weight	28kDa

Antigen Information

Gene Name	RAB7L1
Alternative Names	RAB7L1; Ras-related protein Rab-7L1; Rab-7-like protein 1
Gene ID	8934.0
SwissProt ID	O14966
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB7L1. AA range:115-164


Background

similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitous.,similarity:Belongs to the small


GTPase superfamily. Rab family.,tissue specificity:Ubiquitous.,

Research Area

Image Data



Western blot analysis of lysates from HT-29 cells, using RAB7L1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using RAB7L1 antibody.



Western Blot analysis of various cells using Rab 7L1 Polyclonal Antibody