

Summary

Production Name	Rab L2A Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,IF-P,IF-F,ICC/IF,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RABL2A
Alternative Names	RABL2A; Rab-like protein 2A
Gene ID	11159.0
SwissProt ID	Q9UBK7.The antiserum was produced against synthesized peptide derived from human RABL2A. AA range:179-228

Application

Dilution Ratio	WB 1:500-1:2000, IHC-P 1:100-1:300, ELISA 1:20000, IF-P/IF-F/ICC/IF 1:50-200
Molecular Weight	30-34kDa

Product Name: Rab L2A Rabbit Polyclonal Antibody
Catalog #: AP Rab16764




Background

This gene is a member of the RAB gene family which belongs to the RAS GTPase superfamily. The proteins in the family of RAS-related signaling molecules are small GTP-binding proteins that play important roles in the regulation of exocytotic and endocytotic pathways. This gene maps to the site of an ancestral telomere fusion event and may be a subtelomeric gene. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015], similarity: Belongs to the small GTPase superfamily. Rab family.,

Research Area

Image Data



Western Blot analysis of various cells using Rab L2A Polyclonal Antibody diluted at 1 : 500



Western Blot analysis of COLO205 cells using Rab L2A Polyclonal Antibody diluted at 1 : 500

Note

For research use only.