

Product Name: RAB15 Rabbit Polyclonal Antibody
Catalog #: APRab16775



Summary

Production Name	RAB15 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RAB15
Alternative Names	
Gene ID	376267.0
SwissProt ID	P59190. Synthesized peptide derived from human protein . at AA range: 80-160

Application

Dilution Ratio	WB 1:500-2000, ELISA 1:5000-20000
Molecular Weight	23kDa

Background

function:May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:The GTP bound form of RAB15 interacts


Product Name: RAB15 Rabbit Polyclonal Antibody
Catalog #: AP Rab16775



with REP15.,function:May act in concert with RAB3A in regulating aspects of synaptic vesicle membrane flow within the nerve terminal.,similarity:Belongs to the small GTPase superfamily. Rab family.,subunit:The GTP bound form of RAB15 interacts with REP15.,

Research Area

Image Data



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

Note

For research use only.