

**Product Name: RAB43 Rabbit Polyclonal Antibody****Catalog #: APRab16790**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ELISA
<b>Reactivity</b>	Human,Rat,Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,ELISA 1:5000-1:20000
<b>Molecular Weight</b>	23kDa

**Antigen Information**

<b>Gene Name</b>	RAB43 RAB41
<b>Alternative Names</b>	
<b>Gene ID</b>	339122.0
<b>SwissProt ID</b>	Q86YS6
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 90-170

**Background**

similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.,similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Widely expressed in brain, testis, lung, heart, ovary, colon, kidney, uterus and spleen but not in liver.,

## Research Area

Signaling Pathway; G Protein Signaling; Small G Proteins; Ras Family; Signal Transduction; Protein Trafficking; Nuclear Import / Export; Neuroscience; Processes

## Image Data

