
Product Name: REL2 Rabbit Polyclonal Antibody**Catalog #: APRab17001**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	20kDa

Antigen Information

Gene Name	RLN2
Alternative Names	
Gene ID	6019.0
SwissProt ID	P04090
Immunogen	Synthesized peptide derived from part region of human protein

Background


relaxin 2(RLN2) Homo sapiens This gene encodes a member of the relaxin subfamily and insulin superfamily of peptide hormones. In humans there are three non-allelic relaxin genes. This gene encodes multiple protein isoforms, at least one of which undergoes proteolytic processing. This processing generates relaxin A and B chains that are linked by disulfide bonds to

form the mature peptide hormone. This hormone plays a role in the male and female reproductive systems and was initially noted for its role in pregnancy. This protein also plays broader roles in the cardiovascular system, including in the regulation of blood pressure and control of heart rate, and data from animal models shows that this protein may have anti-fibrotic and cardioprotective effects. [provided by RefSeq, Jul 2016],function:Relaxin is an ovarian hormone that acts with estrogen to produce dilatation of the birth canal in many mammals. May be involved in remodeling of connective tissues during pregnancy, promoting growth of pubic ligaments and ripening of the cervix.,similarity:Belongs to the insulin family.,subunit:Heterodimer of a B chain and an A chain linked by two disulfide bonds.,tissue specificity:Isoform 1 expressed in the ovary during pregnancy. Also expressed in placenta, decidua and prostate. Isoform 2 relatively abundant in placenta, much lower abundance in the prostate gland. Not detected in the ovary.,

Research Area

Signal Transduction

Image Data



Western blot analysis of lysates from HCT116 cells, REL2 Rabbit Polyclonal Antibody was diluted at 1:1000, 4° over night