

Product Name: RNF113B Rabbit Polyclonal Antibody**Catalog #: APRab17284**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:300
Molecular Weight	36kDa

Antigen Information


Gene Name	RNF113B
Alternative Names	RNF113B; RNF161; ZNF183L1; RING finger protein 113B; Zinc finger protein 183-like 1
Gene ID	140432.0
SwissProt ID	Q8IZP6
Immunogen	Synthesized peptide derived from the Internal region of human RNF113B.

Background


similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,

Research Area

Image Data



Western Blot analysis of various cells using RNF113B Polyclonal Antibody



Western Blot analysis of 293 cells using RNF113B Polyclonal Antibody