

Product Name: RNF149 Rabbit Polyclonal Antibody**Catalog #: APRab17289**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	IHC, ICC/IF, ELISA
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000

Molecular Weight

Antigen Information

Gene Name	RNF149
Alternative Names	RNF149; DNAPT2; E3 ubiquitin-protein ligase RNF149; DNA polymerase-transactivated protein 2; RING finger protein 149
Gene ID	284996.0
SwissProt ID	Q8NC42
Immunogen	The antiserum was produced against synthesized peptide derived from human RNF149. AA range: 211-260


Background

domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,domain:The RING-type zinc finger domain mediates binding to an E2 ubiquitin-conjugating enzyme.,pathway:Protein modification; protein ubiquitination.,similarity:Contains 1 PA (protease associated) domain.,similarity:Contains 1 RING-type zinc finger.,

Research Area

Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Ubiquitin E3 Enzymes; RING Finger E3 Ligase

Image Data



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using RNF149 Antibody. The picture on the right is blocked with the synthesized peptide.