
Product Name: RRP22 Rabbit Polyclonal Antibody**Catalog #: APRab17387**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:20000-1:40000
Molecular Weight	30kDa

Antigen Information


Gene Name	RASL10A
Alternative Names	RASL10A; RRP22; Ras-like protein family member 10A; Ras-like protein RRP22; Ras-related protein on chromosome 22
Gene ID	10633.0
SwissProt ID	Q92737
Immunogen	The antiserum was produced against synthesized peptide derived from human RASL10A. AA range:107-156

Background


similarity:Belongs to the small GTPase superfamily. Ras family.,similarity:Belongs to the small GTPase superfamily. Ras family.,

Research Area


Image Data



Western blot analysis of lysates from HeLa cells, using RASL10A Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using RASL10A antibody.



Western Blot analysis of MOUSE-SPLEEN cells using RRP22 Polyclonal Antibody diluted at 1: 1000