

Product Name: SH2D5 Rabbit Polyclonal Antibody**Catalog #: APRab17837**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
Molecular Weight	35kDa

Antigen Information


Gene Name	SH2D5
Alternative Names	SH2D5; SH2 domain-containing protein 5
Gene ID	400745.0
SwissProt ID	Q6ZV89
Immunogen	The antiserum was produced against synthesized peptide derived from human SH2D5. AA range:281-330

Background

similarity:Contains 1 SH2 domain.,similarity:Contains 1 SH2 domain.,

Research Area

Image Data



Western blot analysis of SH2D5 Antibody. The lane on the right is blocked with the SH2D5 peptide.



Western Blot analysis of various cells using SH2D5 Polyclonal Antibody