
Product Name: SLC10A7 Rabbit Polyclonal Antibody**Catalog #: APRab17951**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
Molecular Weight	39kDa

Antigen Information


Gene Name	SLC10A7
Alternative Names	SLC10A7; C4orf13; P7; PSEC0051; Sodium/bile acid cotransporter 7; Na(+)/bile acid cotransporter 7; Solute carrier family 10 member 7
Gene ID	84068.0
SwissProt ID	Q0GE19
Immunogen	The antiserum was produced against synthesized peptide derived from human SLC10A7. AA range:168-217

Background


similarity:Belongs to the sodium:bile acid symporter family.,similarity:Belongs to the sodium:bile acid symporter family.,

Research Area

Image Data



Western blot analysis of lysates from mouse liver, using SLC10A7 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using SLC10A7 antibody.