

---

**Product Name: Slit3 Rabbit Polyclonal Antibody****Catalog #: APRab17974**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	IHC,ICC/IF,ELISA
<b>Reactivity</b>	Human,Mouse,Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

**Dilution Ratio** IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000

**Molecular Weight**

**Antigen Information**

<b>Gene Name</b>	SLIT3
<b>Alternative Names</b>	SLIT3; KIAA0814; MEGF5; SLIL2; Slit homolog 3 protein; Slit-3; Multiple epidermal growth factor-like domains protein 5; Multiple EGF-like domains protein 5
<b>Gene ID</b>	6586.0
<b>SwissProt ID</b>	O75094
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human Slit3.

**Background**


The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration.

Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012],function:May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.,sequence caution:Intron retention.,similarity:Contains 1 CTCK (C-terminal cystine knot-like) domain.,similarity:Contains 1 laminin G-like domain.,similarity:Contains 23 LRR (leucine-rich) repeats.,similarity:Contains 9 EGF-like domains.,tissue specificity:Predominantly expressed in thyroid.,

## Research Area

Axon guidance;

## Image Data



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .