

---

**Product Name: STRN4 Rabbit Polyclonal Antibody****Catalog #: APRab18410**

For research use only.

**Summary**

<b>Description</b>	Rabbit polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ELISA
<b>Reactivity</b>	Human,Mouse
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

**Application**

<b>Dilution Ratio</b>	WB 1:500-1:2000,ELISA 1:5000-1:20000
<b>Molecular Weight</b>	82kDa

**Antigen Information**

<b>Gene Name</b>	STRN4 ZIN
<b>Alternative Names</b>	
<b>Gene ID</b>	29888.0
<b>SwissProt ID</b>	Q9NRL3
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 250-330

**Background**


function: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein., miscellaneous: The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane), similarity: Belongs to the WD repeat striatin family., similarity: Contains 7 WD repeats., function: Binds calmodulin in a

calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane),.similarity:Belongs to the WD repeat striatin family,.similarity:Contains 7 WD repeats.,

## Research Area

Neuroscience

## Image Data



Western Blot analysis of HeLa lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000