
Product Name: TCAL2 Rabbit Polyclonal Antibody**Catalog #: APRab18715**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	24kDa

Antigen Information

Gene Name	TCEAL2 My048 MY0876G0
Alternative Names	
Gene ID	140597.0
SwissProt ID	Q9H3H9
Immunogen	Synthesized peptide derived from part region of human protein


Background

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. [provided by RefSeq, Jul 2008],function:May

be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,

Research Area

Image Data



Western blot analysis of lysates from K562 cells, TCAL2 Rabbit Polyclonal Antibody was diluted at 1:1000, 4° over night