

Product Name: TCEAL7 Rabbit Polyclonal Antibody**Catalog #: APRab18727**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,ELISA
Reactivity	Human,Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
Molecular Weight	

Antigen Information

Gene Name	TCEAL7
Alternative Names	Transcription elongation factor A protein-like 7 (TCEA-like protein 7;Transcription elongation factor S-II protein-like 7)
Gene ID	56849.0
SwissProt ID	Q9BRU2
Immunogen	Synthesized peptide derived from human TCEAL7 AA range: 34-84


Background

function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,function:May be

involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human brain tumor. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .