

Product Name: USP42 Rabbit Polyclonal Antibody**Catalog #: APRab19680**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
Purification	Affinity purification

Application

Dilution Ratio	WB 1:500-1:2000,ELISA 1:5000-1:20000
Molecular Weight	150kDa

Antigen Information

Gene Name	USP42
Alternative Names	USP42; Ubiquitin carboxyl-terminal hydrolase 42; Deubiquitinating enzyme 42; Ubiquitin thioesterase 42; Ubiquitin-specific-processing protease 42
Gene ID	84132.0
SwissProt ID	Q9H9J4
Immunogen	The antiserum was produced against synthesized peptide derived from human USP42. AA range:251-300

Background

catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Broadly expressed.,catalytic activity:Ubiquitin C-terminal thioester + H(2)O = ubiquitin + a thiol.,similarity:Belongs to the peptidase C19 family.,tissue specificity:Broadly expressed.,

Research Area

Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Deubiquitination; Epigenetics and Nuclear Signaling; Ubiquitin & Ubiquitin Like Modifiers

Image Data

