

Product Name: ZNF364 Rabbit Monoclonal antibody**Catalog #: AMRe01416**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Concentration	0.43mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 34 kDa

Antigen Information

Gene Name	RNF115
Alternative Names	E3 ubiquitin-protein ligase RNF115; Rabring 7; Zinc finger protein 364
Gene ID	27246
SwissProt ID	Q9Y4L5
Immunogen	A synthetic peptide of human ZNF364


Background

E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZNF364/RFN115 in PC-12 lysates using ZNF364 antibody.



Western blot analysis of ZNF364/RFN115 in Jurkat, HeLa lysates using ZNF364/RFN115 antibody.