

Summary

Production Name	Rab8A Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	RAB8A
Alternative Names	RAB8A; member RAS oncogene family; MEL; RAB8
Gene ID	4218
SwissProt ID	P61006.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 24 kDa; Observed MW: 24 kDa

Background

Product Name: Rab8A Rabbit Monoclonal Antibody
Catalog #: AMRe01423




The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes.


Research Area

Signal Transduction


Image Data



Western blot analysis of RAB8A in K562, C6, 3T3, HeLa lysates using Rab8A antibody.



Immunocytochemistry analysis of Rab8A(green) in HeLa using Rab8A antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using RAB8A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: Rab8A Rabbit Monoclonal Antibody
Catalog #: AMRe01423



Note

For research use only.