

Product Name: CD32B Rabbit Monoclonal Antibody
Catalog #: AMRe01430



Summary

Production Name	CD32B Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	FCGR2B
Alternative Names	Fc fragment of IgG; low affinity IIb; receptor (CD32); CD32; FCG2; CD32B; FCGR2; IGFR2
Gene ID	2213
SwissProt ID	P31994.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 34 kDa; Observed MW: 44 kDa

Background

Product Name: CD32B Rabbit Monoclonal Antibody
Catalog #: AMRe01430




Receptor for the Fc region of complexed or aggregated immunoglobulins gamma.


Research Area

Immunology

Image Data



Western blot analysis of CD32B in C6 lysates using CD32B antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using CD32B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.