

Product Name: TAF3 Rabbit Monoclonal antibody
Catalog #: AMRe01432



Summary

Production Name	TAF3 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purified

Immunogen

Gene Name	TAF3
Alternative Names	Transcription initiation factor TFIID subunit 3
Gene ID	83860
SwissProt ID	Q5VWG9.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 104 kDa; Observed MW: 140 kDa

Background

Product Name: TAF3 Rabbit Monoclonal antibody
Catalog #: AMRe01432




The TFIID basal transcription factor complex plays a major role in the initiation of RNA polymerase II (Pol II)-dependent transcription.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TAF3 in Jurkat lysates using TAF3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TAF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.