
Product Name: ARFRP1 Rabbit Monoclonal Antibody**Catalog #: AMRe01438**

For research use only.

Summary

Description	Recombinant rabbit monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	0.5mg/ml. The concentration of this product may be batch-dependent.
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:20-1:50
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Antigen Information

Gene Name	ARFRP1
Alternative Names	ARP; Arp1; ARL18
Gene ID	10139
SwissProt ID	Q13795
Immunogen	A synthetic peptide of human ARFRP1


Background

Trans-Golgi-associated GTPase that regulates protein sorting.

Research Area

Signal Transduction

Image Data



Western blot analysis of ARFRP1/ARP in C6 lysates using ARFRP1 antibody.