

Product Name: OXGR1 Rabbit Monoclonal Antibody
Catalog #: AMRe01439



Summary

Production Name	OXGR1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
Purification	Affinity Purification

Immunogen

Gene Name	OXGR1
Alternative Names	P2Y15; Alpha-ketoglutarate receptor 1; G-protein coupled receptor 80; G-protein coupled receptor 99; P2Y purinoceptor 15
Gene ID	27199
SwissProt ID	Q96P68.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 38 kDa; Observed MW: 38 kDa

Product Name: OXGR1 Rabbit Monoclonal Antibody
Catalog #: AMRe01439




Background

Receptor for alpha-ketoglutarate.


Research Area

Signal Transduction

Image Data



Western blot analysis of OXGR1/GPR99 in rat colon lysates using OXGR1 antibody.



Western blot analysis of OXGR1/GPR99 in HeLa lysates using OXGR1/GPR99 antibody.

Note

For research use only.